

# **I. PREANALYTICAL EVENTS**

Clinical Laboratory Medicine  
078 - PEDIATRIC CYTOGENETICS LABORATORY

|  | TECHNICAL<br>CHARGE |
|--|---------------------|
| 078 0005 Amniotic Fluid or Tissue Culture and Chromosomal Analysis | \$788.00            |
| 078 0002 Blood Chromosomal Analysis                                | \$451.00            |
| 078 0006 Cell Culture (Amniotic Fluid or Tissue)                   | \$384.00            |
| 078 0007 Bone Marrow Chromosomal Analysis                          | \$772.00            |
| 078 0008 Blood Chromosomal Analysis of Additional Family Member    | \$385.00            |

Chromosome Anyalysis for Breakage studies-Score 100 cells

Chromosome Analysis; Count additional 50 cells for Mosaicism

Chromosome Analysis for Fragile X

Chromosome Analysis, additional Karyotyping

Specialized Banding Techniques (C, NOR, Q etc.)

In Situ Hybridization for Microdeletions

Chromosome Analysis, Additional High Resolution

## AMNIOTIC FLUID FOR PRENATAL DIAGNOSIS

1. Schedule sample with the Pediatric Cytogenetics Laboratory (409/772-3465).

2. 20-30 cc of amniotic fluid (obtained at 16 weeks post-LMP or 14 weeks post-conception) should be drawn aseptically into a syringe after discarding the first 1-2 cc of fluid. The sample can be delivered in the same syringe capped tightly with a sterile cap, or the sample can be transferred to a sterile tube and capped tightly. Label the sample with the patient's name and date of birth.

3. Complete the Amniotic Fluid Study Request Form (see attached). These will be available in the Private OB Clinic and the Non-Stress Room. Label 3 consult sheets and a Pediatric Cytogenetics Laboratory charge slip.

4. a. UTMB Patient: Sample should be submitted to the laboratory as soon as possible at room temperature packaged well to prevent breakage or leakage.  
b. Off-Campus Patient: Delivery must be within 24 hours of the tap or growth of the cells may be seriously jeopardized.

5. Samples should be delivered to the Pediatric Cytogenetics Laboratory, Suite 3.350, Children's Hospital, in the University of Texas Medical Branch Complex, Galveston, Texas 77555-0359.

6. Preliminary results will be available in 14-18 days. A final written report will be sent in 3-4 weeks. *Physician will be notified by day 7 if there are anticipated delays.*

7. Confirmation of an abnormal amniotic fluid sample result is confirmed through concurrent fetal PUBS blood or subsequent post-mortem skin fibroblast chromosome testing.

Any questions should be directed to either Jerome McCombs, Ph.D., or Lillian H. Lockhart, M.D. at 409/772-3466 or 409/772-3465.

We are in a transitional phase. We have changed our charge slip to reflect the services that we are providing.

The following are examples of how the charge slips should be completed:

1. Stamp or generate charge sticker. This will contain name, UH #, DOB and billing information
2. Check test(s) to be ordered
3. Requesting Faculty
4. Collection Time and Date
5. Receipt time and Date to be added by receiving technologist
6. **It is the responsibility of the technologist receiving the sample to verify the information for completeness and that the paperwork matches the sample.**
7. **Discrepancies are to be recorded in the exceptional case records and reported to the Laboratory Director as soon as possible.**
8. **Inappropriate or incomplete paperwork could be grounds for refusing a sample.**
9. The original copy is retained in the file.
10. The yellow copy is used for laboratory billing and then forwarded to the billing office for professional reimbursement.
11. The new bill provides sample collection and delivery information.
12. Oral request for studies must be followed w/ the appropriate consult sheet and clinical information. Information should be received within 30 days.

# CYTOGENETICS REQUEST FORM

REQUEST FORM TO BE ACCOMPANIED BY TWO STAMPED CONSULT SHEETS.

BLOOD CHROMOSOME ANALYSIS  
5-10MLS OF SODIUM HEPARINIZED BLOOD

AMNIOTIC FLUID OR TISSUE CULTURE AND CHROMOSOME ANALYSIS  
15-30CC AMNIOTIC FLUID - BY PRIOR ARRANGEMENT WITH LAB

TISSUE CULTURE AMNIOTIC FLUID OR TISSUE CULTURE  
BY PRIOR ARRANGEMENT WITH LAB

BONE MARROW CHROMOSOME ANALYSIS  
BY PRIOR ARRANGEMENT WITH LAB

BLOOD CHROMOSOME ANALYSIS; ADDITIONAL FAMILY MEMBER  
BY PRIOR ARRANGEMENT WITH LAB

TESTS BELOW MUST INCLUDE ONE OF THE ABOVE SERVICES

CHROMOSOME ANALYSIS FOR BREAKAGE STUDIES - SCORE 100 CELLS

CHROMOSOME ANALYSIS; COUNT ADDITIONAL 50 CELLS FOR MOSAICISM

CHROMOSOME ANALYSIS FOR FRAGILE X

CHROMOSOME ANALYSIS, ADDITIONAL KARYOTYPING

SPECIALIZED BANDING TECHNIQUES (C,NOR,Q etc.)

IN SITU HYBRIDIZATION FOR MICRODELETIONS

CHROMOSOME ANALYSIS, ADDITIONAL HIGH RESOLUTION

DELIVER TO: CYTOGENETICS LABORATORY  
3.350 CHILDRENS HOSPITAL  
EXTENSIONS: 23466 OR 23465



|   |                           |   |
|---|---------------------------|---|
| REQUESTING<br>PHYSICIAN:  | DIAGNOSIS:<br>ICD-9 CODE: | COLLECTED BY:<br>DATE / TIME:   |
| Unit History #:<br>Name:<br>Financial #:<br>Location:   |                           | Results to Physician (Name and ID #):<br><br>Floor / Clinic Phone Extension:<br><br>UTMB Doctor # / Beeper #: |
| <b>UNIVERSITY OF TEXAS MEDICAL BRANCH</b><br>LABORATORY MEDICINE - GALVESTON, TX 77555-0743<br>CLIA #45D0660281 CAP #21251-01<br>MICHAEL BISSEL, M.D., Ph.d., M.P.H. - CLINICAL DIRECTOR<br>ABIDA K. HAQUE, M.D. - ANATOMIC DIRECTOR<br>TELEPHONE: (800) LAB 2266 |                           |   |

PLACE INVISION LABEL HERE

PLACE LIS LABEL HERE

CY 0001202

|         |                 |         |                 |
|---------|-----------------|---------|-----------------|
| NAME:   | NAME:           | NAME:   | NAME:           |
| UH#: CY | UH#: CY 0001202 | UH#: CY | UH#: CY 0001202 |

# CYTOGENETICS REQUEST FORM

REQUEST FORM TO BE ACCOMPANIED BY TWO STAMPED CONSULT SHEETS.

BLOOD CHROMOSOME ANALYSIS  
5-10MLS OF SODIUM HEPARINIZED BLOOD

AMNIOTIC FLUID OR TISSUE CULTURE AND CHROMOSOME ANALYSIS  
15-30CC AMNIOTIC FLUID - BY PRIOR ARRANGEMENT WITH LAB

TISSUE CULTURE AMNIOTIC FLUID OR TISSUE CULTURE  
BY PRIOR ARRANGEMENT WITH LAB

BONE MARROW CHROMOSOME ANALYSIS  
BY PRIOR ARRANGEMENT WITH LAB

BLOOD CHROMOSOME ANALYSIS; ADDITIONAL FAMILY MEMBER  
BY PRIOR ARRANGEMENT WITH LAB

TESTS BELOW MUST INCLUDE ONE OF THE ABOVE SERVICES

CHROMOSOME ANYALYSIS FOR BREAKAGE STUDIES - SCORE 100 CELLS

CHROMOSOME ANALYSIS; COUNT ADDITIONAL 50 CELLS FOR MOSAICISM

CHROMOSOME ANALYSIS FOR FRAGILE X

CHROMOSOME ANALYSIS, ADDITIONAL KARYOTYPING

SPECIALIZED BANDING TECHNIQUES (C,NOR,Q etc.)

IN SITU HYBRIDIZATION FOR MICRODELETIONS

CHROMOSOME ANALYSIS, ADDITIONAL HIGH RESOLUTION

DELIVER TO: CYTOGENETICS LABORATORY  
3.350 CHILDRENS HOSPITAL  
EXTENSIONS: 23466 OR 23465

REQUESTING

PHYSICIAN:

DIAGNOSIS:

ICD-9 CODE:

COLLECTED BY:

DATE / TIME:

Unit History #:

Name:

Financial #:

Location:

UNIVERSITY OF TEXAS MEDICAL BRANCH  
LABORATORY MEDICINE - GALVESTON, TX 77555-0743  
CLIA #45D0660281 CAP #21251-01  
MICHAEL BISSEL, M.D., Ph.d., M.P.H. - CLINICAL DIRECTOR  
ABIDA K. HAQUE, M.D. - ANATOMIC DIRECTOR  
TELEPHONE: (800) LAB 2266

Results to Physician (Name and ID #):

Floor / Clinic Phone Extension:

UTMB Doctor # / Beeper #:

PLACE INVISION LABEL HERE

PLACE LIS LABEL HERE

CY0001202

**SAMPLE PRIORITY** (Most to Least Crucial)

1. Percutaneous Blood Sampling (PUBS)\*
2. Peripheral blood from newborn
3. Peripheral blood from pregnant couple
4. Amniotic fluid
  - a. 20 week's gestation
  - b. Less than 20 weeks
  - c. More than 20 weeks
5. Bone marrow
6. Other peripheral blood samples
  - a. Non-fragile(x), non-Fanconi's anemia
  - b. Fragile(x)
  - c. Fanconi's anemia
7. Skin fibroblasts
  - a. Metabolic studies
  - b. Autopsy specimens
8. Research samples

\*PUBS are to be reported (preliminary) within 72 hrs with appropriate documentation in the chart. The written report is to be finalized within seven days.

\* Any delay that can be anticipated that would

## AMNIOTIC FLUID FOR PRENATAL DIAGNOSIS

1. Schedule sample with the Pediatric Cytogenetics Laboratory (409/772-3465).

2. 20-30 cc of amniotic fluid (obtained at 16 weeks post-LMP or 14 weeks post-conception) should be drawn aseptically into a syringe after discarding the first 1-2 cc of fluid. The sample can be delivered in the same syringe capped tightly with a sterile cap, or the sample can be transferred to a sterile tube and capped tightly. Label the sample with the patient's name and date of birth.

3. Complete the Amniotic Fluid Study Request Form (see attached). These will be available in the Private OB Clinic and the Non-Stress Room. Label 3 consult sheets and a Pediatric Cytogenetics Laboratory charge slip.

4. a. UTMB Patient: Sample should be submitted to the laboratory as soon as possible at room temperature packaged well to prevent breakage or leakage.  
b. Off-Campus Patient: Delivery must be within 24 hours of the tap or growth of the cells may be seriously jeopardized.

5. Samples should be delivered to the Pediatric Cytogenetics Laboratory, Suite 3.350, Children's Hospital, in the University of Texas Medical Branch Complex, Galveston, Texas 77555-0359.

6. Preliminary results will be available in 14-18 days. A final written report will be sent in 3-4 weeks. *Physician will be notified by day 7 if there are anticipated delays.*

7. Confirmation of an abnormal amniotic fluid sample result is confirmed through concurrent fetal PUBS blood or subsequent post-mortem skin fibroblast chromosome testing.

Any questions should be directed to either Jerome McCombs, Ph.D., or Lillian H. Lockhart, M.D. at 409/772-3466 or 409/772-3465.



PEDIATRIC CYTOGENETICS LABORATORY  
331 CHILDREN'S HOSPITAL, UTMB, GALVESTON, TEXAS 77555-0359  
409/772-3465

**AMNIOTIC FLUID STUDY REQUEST FORM**

**PLEASE DISCARD FIRST 1-2 CC OF AMNIOTIC FLUID**

Date \_\_\_\_\_ Time \_\_\_\_\_ (Lab Use Only - Lab No \_\_\_\_\_ )

Patient's Name \_\_\_\_\_ UH#: \_\_\_\_\_

Patient's Address \_\_\_\_\_

Patient's Phone \_\_\_\_\_ Date of Birth \_\_\_\_\_ LNMP \_\_\_\_\_

Gestational Age: By ultrasound \_\_\_\_\_ When \_\_\_\_\_ By Dates \_\_\_\_\_

Name of UTMB Physician Performing Amniocentesis \_\_\_\_\_

What is the reason for the amniocentesis? \_\_\_\_\_

NAME AND ADDRESS OF OB SATELLITE CLINIC PATIENT ATTENDS \_\_\_\_\_

IF PRIVATE, NAME AND ADDRESS OF PRIVATE OBSTETRICIAN PATIENT  
ATTENDS \_\_\_\_\_

\_\_\_\_\_ Telephone \_\_\_\_\_

DID YOU DISCARD THE FIRST 1-2 CC OF AMNIOTIC FLUID? Yes \_\_\_\_\_ No \_\_\_\_\_

If not, why \_\_\_\_\_

**LABORATORY USE ONLY**

Number of syringes received \_\_\_\_\_

Color of fluid \_\_\_\_\_ Contaminated with blood? Yes \_\_\_\_\_ No \_\_\_\_\_

Collected \_\_\_\_\_ Received \_\_\_\_\_ Culture initiated \_\_\_\_\_

Cell pellet description (RBC's, size of pellet, etc.) \_\_\_\_\_

Number of flasks/coverslips set up \_\_\_\_\_ Tissue culture medium/  
antibiotics used \_\_\_\_\_

**THIS FORM MUST BE COMPLETED AND SENT WITH EACH AMNIOTIC  
FLUID SPECIMEN RECEIVED BY THE PEDIATRIC CYTOGENETICS  
LABORATORY (CHARGE SLIP AND 3 LABELED CONSULTS ALSO)**

## PERIPHERAL BLOOD FOR CHROMOSOME STUDIES

1. Schedule sample with the Pediatric Cytogenetics Laboratory (409/772-3465).
2. Blood must be heparinized and sterile. Either wet the inside of a syringe with sodium heparin (without preservatives) leaving 0.10 cc before drawing the blood (2.0-5.0 ml); mix well; label sample with patient's name and date of birth OR draw blood (2.0-5.0 ml) into a vacutainer containing sodium heparin (without preservatives); mix well, label sample with patient's name and date of birth.
3.
  - a. UTMB Patient: Label 3 consult sheets and one Pediatric Cytogenetics Laboratory charge slip and provide appropriate information.
  - b. Off-Campus Patient: Contact laboratory for more information.
4. Sample should be submitted to the laboratory as soon as possible at room temperature. If there is a delay in delivery (overnight or weekend), the sample may be placed in the refrigerator (but not near the freezer). Keep it cool but not frozen. If delivering the sample during the summer, avoid overheating. Placing the sample in a cup of water with a small amount of ice may be helpful.
5. Samples should be delivered to the Pediatric Cytogenetics Laboratory, Room 331, Children's Hospital, in the University of Texas Medical Branch Complex, Galveston, Texas 77555-0359.
6. Preliminary results will be available in 7-10 days. A final written report will be sent in 3-4 weeks.
7. Any questions should be directed to either Jerome McCombs, Ph.D. or Lillian H. Lockhart, M.D. at 409/772-3466 or 409/772-3465.
8. A copy of the Laboratory Survival Manual is available at all nursing stations. The sampling procedures therein are reviewed annually.

## BONE MARROW FOR CHROMOSOME STUDIES

1. Schedule sample with the Pediatric Cytogenetics Laboratory (409/772-3465).
2. Collect 2-5 ml of aspirate into: Either a vacutainer with sodium heparin OR a sodium heparinized syringe. Label sample with patient's name and date of birth.
3.
  - a. UTMB Patient: Label 3 consult sheets and one Pediatric Cytogenetics Laboratory charge slip and provide appropriate information.
  - b. Off-Campus Patient: Contact laboratory for more information.
4.
  - a. UTMB Patient: Call laboratory and technician will come to pick up the sample.
  - b. Off-Campus Patient: Sample should be submitted to the laboratory as soon as possible at room temperature, packaged well to prevent breakage or leakage. Delivery must be within 24 hours of collection or quality of the specimen may be greatly diminished. Store in the refrigerator until the sample is shipped. Do not freeze.
5. Samples should be delivered to the Pediatric Cytogenetics Laboratory, Suite 3.350, Children's Hospital, in the University of Texas Medical Branch Complex, Galveston, Texas 77555-0359. We can also arrange to pick up the sample ourselves.
6. Preliminary results will be available in 7-10 days. A final written report will be sent in 3-4 weeks.
7. Any questions should be directed to either Jerome McCombs, Ph.D. or Lillian H. Lockhart, M.D., at 409/772-3466 or 409/772-3465.

### SKIN SAMPLE FOR CHROMOSOME STUDIES

1. Schedule sample with the Pediatric Cytogenetics Laboratory (409/772-3465).
2. A roughly 4 mm<sup>2</sup> section of skin is taken from the deceased infant by the Cytogenetics technician using sterile techniques.
3. A consult sheet and Pediatric Cytogenetics charge slip are filled out by Autopsy Service.
4. The sample is returned to the Pediatric Cytogenetics Laboratory, Suite 3.350, Children's Hospital, UTMB, Galveston, Texas 77555-0359 and put to culture immediately.
5. Preliminary results will be available in roughly four weeks. A final written report will be sent in 5-6 weeks.
6. Any questions should be directed to either Jerome McCombs, Ph.D., or Lillian H. Lockhart, M.D., at 409/772-3466.

### SKIN SAMPLE FOR METABOLIC STUDY

1. Schedule sample with the Pediatric Cytogenetics Laboratory (409/772-3465). Pick up a sterilized skin pack and test tube with growth medium from the Cytogenetics lab.
2. Using the skin punch in the skin pack, collect a skin sample. Label sample with patient's name and date of birth.
3. Submit sample to the Pediatric Cytogenetics Laboratory and return skin pack with consult sheets and Pediatric Cytogenetics Lab and Developmental Nutrition and Metabolism Laboratory charge slips. Deliver to Suite 3.350, Children's Hospital, UTMB, Galveston, Texas 77555-0359.
4. Preliminary results will be available from Developmental Nutrition and Metabolism Laboratory in 6-8 weeks. A final written report will follow shortly.
5. Any questions should be directed to either Jerome McCombs, Ph.D., 409/772-3466 or David Rassin, Ph.D. 409/772-1139.

## INITIAL SAMPLE ACCEPTANCE AND DATA REQUIREMENTS

All samples accepted must have the following information:

1. Labeled sample tube with Patient's name, UH#, or corresponding labeled charge slip number.
2. Three labeled consult sheets with Physician's name, phone number, and reason for referral.
3. Samples will be accepted without consult sheets only if the referring physician is identified and a verbal reason for referral is given.

Time and date the sample is received is printed on the charge slip as illustrated on the following page.

Samples not meeting minimal standards. These are to be called to the laboratory director's attention immediately. At that time a decision will be made by the director as to acceptance/rejection and possible corrective actions are listed below: **(Note: Universal safety precautions in a biological safety cabinet are employed when setting up samples.)**


1. **Amniotic Fluid**
  - A. Bloody--when growth established, flask washed and new medium added.
  - B. <15 ml received--See sheet.
  - C. If doctor did not discard first few cc's, watch for maternal cell contamination. Contact physician and inform him of potential diagnostic difficulties.
  - D. HIV+, hepatitis+, venereal disease+--double gloves, discard waste in separate biohazard bag, pour fluid waste in container with bleach.
2. **Blood**
  - A. Blood in contact with lithium heparin or EDTA--spin out cells and resuspend in medium with sodium heparin right away.

- B. Blood syringe on ice--remove from ice immediately.
- C. Blood received on Wednesdays:
  - 1. Set up Thursday for 96-hour culture
  - 2. Can set up 72-hour synchronized culture on Friday.
    - a. Store syringe in refrigerator
- D. HIV+, hepatitis+, venereal disease+--see Amniotic Fluid step D above.

Any sample problems or unusual conditions are noted in the patient's chart and reported to laboratory director immediately.
- 3. **Bone Marrow**
  - A. <1 ml received--still set up 2 cultures (direct and overnight)
  - B. No visible spicules--proceed but may get less desirable results
  - C. HIV+, hepatitis+, venereal disease+--see Amniotic Fluid step D above
  - D. BM received on Fridays--do directs
- 4. **Skin** (morgue specimens, specifically)
  - A. If baby in formalin, necrotic skin, dead more than a day: take organ tissue (e.g., heart, lung)
  - B. HIV+, hepatitis+, venereal disease+ (baby or its mother): don mask and apron in addition to gloves when taking sample--Refer to Amniotic Fluid step D above

**Note:** A sheet outlining the unacceptable status of a sample is placed in the patient's chart.

Approved: \_\_\_\_\_



Sample Collection

---

### **PROTOCOL FOR ACCEPTING < 15 ML AMNIOTIC FLUID**

1. Inform lab director that < 15 ml of amniotic fluid was received.
2. Contact referring physician to identify any possible problem that could affect the cell culture and to confirm that the physician is aware of our amniotic fluid requirement of 20-33 cc.
3. Initiate two flasks (25 cm).
4. Perform standard analysis of 20 cells.
5. If the standard analysis reveals a 46,XX, count an additional 50 cells with fluorescence (Q-bands).
6. State in the report (if female is identified) that due to the limited sample the possibility exists for maternal cell contamination and this cannot be excluded.



## GUIDELINES FOR LABELING PATIENT SAMPLE TUBES, SLIDES, CAPTURED IMAGES, AND PHOTOGRAPHS

1. Samples are initiated following standard procedures in the protocol manual.
2. Only one patient is initiated (cultures started in biosafety cabinet) at a time.
3. Place on each flask a label (tape) containing:
  - A. Name
  - B. Laboratory #
  - C. Date initiated
  - D. Type medium
  - E. + or - EB
4. When transferring sample to conical centrifuge tube, merely peel tape from flask and place on centrifuge tube.
5. Color code pipets that are used with each sample to match sample tape.
6. When slides are made, each slide is to be labeled as it is made. Label will include:
  - A. Name
  - B. Laboratory #
  - C. Type medium
  - D. + or - EB/Synch
7. Slides are to be prepared on one patient at a time.
8. It has been out policy to retain all slides, negatives, and prints permanently on all cases, normal and abnormal. Note: In APRIL, 1996, a PSI imaging system completely replaced our darkroom facility. CAPTURED IMAGES have legends containing name, number, date and slide # and coordinates. All karyotypes generated continue to show patient name, number, karyotype designation, and karyotype date.

Dr. \_\_\_\_\_  
Address: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

☐ IN ADDITION TO WRITTEN REPORT, PLEASE CALL EXT / BEEPER# BELOW.

TO: PHYSICIAN OR SERVICE

Reason for Requesting Consultation – This request is: ☐ URGENT ☐ ROUTINE

Signature:  
(Physician)

Ext.#

Beeper:

Date:

REPORT: \_\_\_\_\_

I have examined the patient and reviewed the medical record.

Signature/Title: \_\_\_\_\_ Ext. # \_\_\_\_\_ Beeper: \_\_\_\_\_

Service: \_\_\_\_\_ Date: \_\_\_\_\_

## DATA ENTRY AND CODING

All samples are entered sequentially into the Patient Log Book located in the Microscope Room.

### Data to be entered:

Lab #[sequential # followed by suffix indicating sample type (BL-blood; AF- amniotic fluid; SK-skin; BM-bone marrow)]

Date Received

Patient Name

UH Number

DOB

Race

Service

Pay Class

Charge Slip Sent

Referring Physician

Indication

Date Completed

ISCN Diagnosis

Additional Information (e.g., AF, AFP results, etc.)

Patient/sample information entered on daily tissue culture log sheets (ring binder in microscope room), AF log sheets filed in patient's folder when study is completed.

The following table lists indication for cytogenetic studies and possible chromosomal causes. This is merely a general guideline and does not preclude a complete analysis.

## REASON FOR REFERRAL

| <u>Problem</u>                | <u>Possible Anomalies</u>   |
|-------------------------------|-----------------------------|
| Ambiguous Genitalia           | X and Y - ACCBI             |
| Amenorrhea                    | X                           |
| Aniridia                      | 11p13                       |
| Azoospermia                   | X and Y - ACCBI             |
| Behavior Problems             | X and Y - ACCBI             |
| Cat-Eye Syndrome              | 22pter-q11                  |
| Cleft Lip and Palate          | ACCBI                       |
| Congenital Heart Disease      | ACCBI                       |
| Cornelia-de-Lange             | tri 3q or 13q               |
| Developmental Delay           | ACCBI                       |
| Di George                     | 22q11                       |
| Epilepsy                      | ACCBI                       |
| Exomphalos                    | 18-ACCBI                    |
| Failure to Thrive             | ACCBI                       |
| Floppy Baby                   | 15q12,+21 - ACCBI           |
| Gynecomastia                  | X or Y                      |
| Hypospadia                    | del of ACCBI                |
| Infertility                   | X or Y, Translocation-ACCBI |
| Langer-Giedions Syndrome      | 8q22-24                     |
| Macro-orchidism               | fra(X)                      |
| Mental Retardation            | fra(X)-ACCBI                |
| Microcephaly                  | 4p, 5p, 13q-,18-ACCBI       |
| Multiple Congenital Anomalies | ACCBI                       |
| Odd Facies                    | ACCBI                       |
| Prader-Willi                  | 15q11-12                    |
| Recurrent Miscarriages        | ACCBI                       |
| Retinoblastoma                | 13q14                       |
| Rocker Bottom Feet            | 13 or 18                    |
| Small Stature                 | X                           |
| Seizures                      | 4p, 9p-, ACCBI              |
| Tall                          | X or Y                      |
| Wilms Tumor                   | 11p13                       |

\* Any Chromosome Could Be Involved